

常州健一

Changzhou Health Import and Export Company Ltd

常州健一进出口有限公司

https://www.healthtopmed.com/







COMPANY PROFILE

Founded in 2019, CHANGZHOU HEALTH IMPORT AND EXPORT COMPANY LTD, a young and aggressive company mainly engaged in medical product development, covering thousands of medical products, is dedicating to becoming a market leader in the medical products. Our company is known for its good quality reputation, moderate price, and thoughtful services.

We are a professional manufacturer mainly produce various disposable medical products, such as disposable syringes, auto-destroy syringes, insulin syringes, oral syringes, hypodermic needles, infusion & transfusion sets, IV catheter, and all other kinds of medical products.

Based on our principles of honesty and joint efforts among our staff, our company has been constantly expanding to take a leading position in the medical industry. Moreover, our highly efficient group of professional technicians has ceaselessly put forth new products every year, thus maintaining the robust growth trend of the company.

In order to strengthen our management, as well as to ensure that such important products as those within the medical industry meet global quality requirements, we have attained ISO13485 and CE certification for most of our products.We aim to high quality, safety and availability of medical products supplied.

Our principles:

We are working to help people enjoy health.

We are working to use our knowledge and experience to help the patients get well.

We are working to create the best working conditions in which the best employees can develop and achieve their goals.

We are working to strengthen the leadership status of our company in the pharmaceutical market.



Δ

4

Characterian Ch

iso

Recisely Right.

Attachment to Cardinate Registration No.: DD 60144633 0001 Report No.: 15092120 002 Manufacturer: Changahan Yangia Medical Appliance Co., LSL Bandhatos Denet, Zhengja Ton Taching Distoit Taching Distoit

 Paport No.
 2405402011-0

 Effective cales.
 2025-85-06

 Charlow cales.
 2023-81-06

 Move data
 2023-81-06

 Move data
 2023-81-06

Application for QMB Devicements SA 2014718

50

CE

CONTENTS 目录



Syringes & Needles

Disposable Auto-destroy Syringe	1
Disposable Sterile Syringe 3-part	
Insulin Syringe	
Disposable Sterile Syringe 2-part	
Enteral Feeding Syringe	4
Irrigation Syringe	4
Oral Syringe	4
Sterile Hypodermic Needle	
Disposable Scalp Vein Needle	5
Blood Safety Lancet	6
Disposable Blood Lancet	6
Blood Collection Needle	6
Spinal Needle	7
Epidural Needle	7

□ Infusion & Transfusion

Infusion Sets ·····	8-9
Transfusion Sets	10
Infusion Sets With Burette	10
Extension Tube	10
Three way Stopcock	10
IV Cannula/Catheter	11

Emergency & Surgery

Tourniquet	
Surgical Blades	



Other Medical Consumables

Alcohol Pad	
Tongue Depressor12	
Cotton Applicator12	

Certificate



health



DISPOSABLE AUTO-DESTROY SYRINGE Capacity: 1ml,2ml,3ml,5ml,10ml

Ref. No.	Needle Size	Remark
TD -AD -A -0.5ml		
TD -AD -A -1ml	18G -30G A	
TD -AD -A -2ml		A: With assembled needle
TD -AD -A -3ml		
TD -AD -A -5ml		
TD -AD -A -10ml		
TD -AD -B -0.5ml		
TD -AD -B -1ml		
TD -AD -B -2ml	18G -30G	B: With needle at side
TD -AD -B -3ml	18G -30G	
TD -AD -B -5ml		
TD -AD -B -10ml		

PASE DAS

DISPOSABLE AUTO-DESTROY SYRINGE Capacity: 0.5ml

Instruction For Use:

1.Pull back the plunger to charge the syringe with solution.

2.Press plunger forward to complete injection until it reaches the stop position. Lock mechanism will be activated and locks plunger in stop position. An audible "Click" will sound.

3. Forcing the plunger backward will cause it to break. Dispose safely in a disposable container



SYRINGES & NEEDLES



DISPOSABLE STERILE SYRINGE 3-PART LUER SLIP

Capacity: 1ml,2ml,2.5ml,3ml,5ml,10ml,20ml,30ml,50ml,60ml

Ref. No.	Needle Size	Remark	Ref. No.	Needle Size	Remark
TD -SN -A -1ml			TD -SN -C -1ml		
TD -SN -A -2ml			TD -SN -C -2ml		
TD -SN -A -2.5ml			TD -SN -C -2.5ml		
TD -SN -A -3ml	18G -30G		TD -SN -C -3ml		
TD -SN -A -5ml		A: Centrical Luer Slip with	TD -SN -C -5ml	18G -30G	C: Centrical
TD -SN -A -10ml		Assembled Needle	TD -SN -C -10ml	186 -306	Luer Slip with Needle at Side
TD -SN -A -20ml			TD -SN -C -20ml		
TD -SN -A -30ml			TD -SN -C -30ml		
TD -SN -A -50ml			TD -SN -C -50ml		
TD -SN -A -60ml			TD -SN -C -60ml		
TD -SN -B -10ml			TD -SN -D -10ml		
TD -SN -B -20ml	18G -30G	B: Excentrical	TD -SN -D -20ml		
TD -SN -B -30ml		Luer Slip with Assembled	TD -SN -D -30ml	18G -30G	D: Excentrical Luer Slip with Needle at Side
TD -SN -B -50ml		Needle	TD -SN -D -50ml		ut oldo
TD -SN -B -60ml			TD -SN -D -60ml		



DISPOSABLE STERILE SYRINGE 3-PART LUER LOCK

Capacity: 1ml,2ml,2.5ml,3ml,5ml,10ml,20ml,30ml,50ml,60ml

Blister pack / Poly bag With/Without needle,With/ Without cap Natural latex / Synthetic rubber(latex free)

Ref. No.	Needle Size	Remark	Ref. No.	Needle Size	Remark
TD -SN -E -1ml			TD -SN -F -1ml		
TD -SN -E -2ml			TD -SN -F -2ml		
TD -SN -E -2.5ml			TD -SN -F -2.5ml		
TD -SN -E -3ml			TD -SN -F -3ml		
TD -SN -E -5ml	18G -30G	E: Luer Lock with Assembled Needle	TD -SN -F -5ml	18G -30G	F: Luer Lock
TD -SN -E -10ml	186 -306		TD -SN -F -10ml	186 -306	with Needle at Side
TD -SN -E -20ml			TD -SN -F -20ml		
TD -SN -E -30ml			TD -SN -F -30ml		
TD -SN -E-50ml			TD -SN -F -50ml		
TD -SN -E-60ml			TD -SN -F -60ml		



SYRINGES & NEEDLES



INSULIN SYRINGE Capacity: 0.3ml,0.5ml,1ml

Ref. No.	Needle Size	Remark
TD -INS -A -0.3ml -U100		
TD -INS -A -0.5ml -U100		
TD -INS -A -1ml -U100	29G -31G	 A: fixed needle, with protective end cap
TD -INS -A -0.5ml -U40		
TD -INS -A -1ml -U40		
TD -INS -B -0.3ml -U100		
TD -INS -B -0.5ml -U100		
TD -INS -B -1ml -U100	29G -31G	B:fixed needle, without protective end cap
TD -INS -B -0.5ml -U40		
TD -INS -B -1ml -U40		
TD -INS -C -1ml -U100	26G -31G	C: with a 6 % (Luer) male conical fitting and a detachable assembled needle
TD -INS -D -1ml -U100	26G -31G	D: with a 6 % (Luer) male conical fitting and a detachable needle at side

Sterile, Blister pack / Poly bag

Ref. No.	Needle Size	Remark
TD -STP -A -2ml		A:Centrical Luer
TD -STP -A -5ml		Slip
TD -STP -B -5ml	18G -30G	
TD -STP -B -10ml		B: Excentrical Luer Slip
TD -STP -B -20ml		



DISPOSABLE STERILE SYRINGE 2-PART

Capacity: 2ml,5ml,10ml,20ml



SYRINGES & NEEDLES





Product Name	Ref. No.	
Enteral Feeding Syringe	TD -EFS -20ml	
	TD -EFS -30ml	
	TD -EFS -50ml	
	TD -EFS -60ml	

ENTERAL FEEDING SYRINGE Capacity: 20ml,30ml,50ml,60ml

Product Name	Ref. No.
Irrigation Syringe	TD -CS -20ml
	TD -CS -50ml
	TD -CS -60ml
	TD -CS -100ml

IRRIGATION SYRINGE

Capacity: 20ml,50ml,60ml,100ml



Product Name	Ref. No.
Oral Syringe	TD -OS -1ml
	TD -OS -3ml
	TD -OS -5ml
	TD -OS -10ml

ORAL SYRINGE

Capacity: 1ml,3ml,5ml,10ml

INSTRUCTIONS FOR USE

For single patient use only

This medicine syringe is suitable for oral, liquid dose medication only

1.Fit the bottle adaptor firmly into the neck of medicine bottle 2.Empty syringe and draw out plunger to dosage mark required.

3.Fit syringe into bottle adaptor and invert bottle.

4.Push in plunger fully and then slowly draw out medicine to dosage mark required.

5.Check for any air bubbles in the syringe, if any are present repeat step 4 until bubbles disappear.

6.To measure dose accurately, line the end back ring, on plunger with the required dosage mark.

7.Place medicine bottle upright and remove syringe. Check dosage accurately again.

8.Check patient is sitting or held upright before giving medicine. 9.Place syringe towards inside of cheek and slowly push the plunger, giving the patient time to swallow.

Rapid squirting of medicine can cause choking.

10.Remove bottle adaptor and replace cap securely on medicine bottle.







health

STERILE HYPODERMIC NEEDLE

Si	ze	Product	Co	lor	Length					
OD(mm)	Gauge	Code	Color	&Code	1 / 2''	5/8''	3 / 4''	1''	1 1/4''	1 1/2"
0.30	30G	HN30G		light yellow						
0.33	29G	HN29G		red						
0.36	28G	HN28G		blue green						
0.40	27G	HN27G		grey						
0.45	26G	HN26G		brown						
0.50	25G	HN25G		orange						
0.55	24G	HN24G		purple						
0.60	23G	HN23G		blue						
0.70	22G	HN22G		black						
0.80	21G	HN21G		green						
0.90	20G	HN20G		yellow						
1.10	19G	HN19G		beige						
1.20	18G	HN18G		pink						
1.60	16G	HN16G		white						
1.80	15G	HN15G		blue grey						

Technical Specification:

Needle Gauge: 15G-30G

Needle Material: Stainless Steel

Needle Hub Material: PP

Sterilization: Sterilized by EO gas

Shelf Life: 5 years






Technical Specification:

Luer slip / Luer lock Needle Gauge: 18G- 27G Needle Material: Stainless Steel Tube Material: PVC Sterilization: Sterilized by EO gas Shelf Life: 5 years

DISPOSABLE SCALP VEIN NEEDLE



SYRINGES & NEEDLES





BLOOD SAFETY LANCET Needle available:21G,23G,26G,28G,30G

Gamma Radiation



DISPOSABLE BLOOD LANCET

Needle available:21G,23G,26G,28G,30G,33G Gamma Radiation

BLOOD COLLECTION NEEDLE

Needle available:16G,18G,20G,21G,22G,23G





health

OLINCKE	POINT TYPE	PENCI	POINT TYPE
QUINCILL			

SPINAL NEEDLE

Product		Length	
Product	mm	Gauge	mm
Quincke Point Type	0.4-1.6	27G-16G	25-200

Completely transparent needles stand provides convenience to observe the outflow of the cerebrospinal fluid.

The needle core with pin allows a good joint of needle edge.

Product		Length	
Product	mm	Gauge	mm
Pencil Point Type	0.4-0.9	27G-20G	25-200

Specially designed needlepoint shape causes no trauma of spinal dura mater, and automatically closes the puncture port after withdrawal of the needle to reduce the leak of cerebrospinal fluid, with the occurrence of PDPH evidently declined.

Proper shape and position of the side hole maintains adequate flow and reduces imperfect retardations as well.

A		

Product		Length	
Product	mm	Gauge	mm
Epidural Needle	0.9 -1.8	20G -15G	40 -150

Bend and round head helps to reduce largely the risk of piercing spinal dura mater and to ensure the smooth pass of the canula.

The needle tube is scaled for the purpose of depth control.

The needle stand with wings is easy to hold and to perform needle insertion.

EPIDURAL NEEDLE



INFUSION & TRANSFUSION



Ref. No.	Specification
TD -ISG -A	Closure -piercing device, integral air - inlet with air filter and closure, drip chamber, flow regulator, latex bulb, luer slip

Ref. No.	Specification
TD -ISG -B	Closure –piercing device, integral air – inlet with air filter and closure, drip chamber, flow regulator, latex bulb, luer lock



Ref. No.	Specification
TD -ISG -C	Closure -piercing device, integral air - inlet with air filter and closure, drip chamber, flow regulator, Y-port, luer slip



Ref. No.	Specification
TD -ISG -D	Closure -piercing device, integral air - inlet with air filter and closure, drip chamber, flow regulator, Y-port, luer lock



INFUSION & TRANSFUSION



Ref. No.	Specification
TD -ISG -F	Closure -piercing device, integral air - inlet with air filter and closure, drip chamber, flow regulator, luer lock

	A	

Ref. No.	Specification		
TD -ISG -P	Closure -piercing device, integral air - inlet with air filter and closure, drip chamber, flow regulator, Y-port, luer slip with cap		

Specification

inlet with air filter and closure, drip chamber, flow regulator, luer slip


Ref. No.

TD -ISG -E

Ref. No.	Specification		
TD -ISG -N	Closure –piercing device, integral air – inlet with air filter and closure, drip chamber, flow regulator, Y–port, luer lock with cap		



INFUSION & TRANSFUSION





INFUSION & TRANSFUSION





IV CANNULA INJECTION TYPE Size: 14G,16G,18G,20G,22G,24G

IV CANNULA WING TYPE Size: 14G,16G,18G,20G,22G,24G

Si	ze	Color	Length					
OD(mm)	Gauge	Color	3/4''	1''	1 1/2''	1 1/4''	2	2 1/2''
2.00	14G	Orange						
1.70	16G	Grey						
1.20	18G	Green						
1.00	20G	Pink						
0.90	22G	Blue						
0.70	24G	Yellow						
0.60	26G	Purple						



EMERGENCY & SURGERY





TOURNIQUET

Size: 2.5*45cm

Color: Dark / Light blue, orange, green, yellow, pink and more




SURGICAL BLADES Stainless steel / Carbon steel Gamma radiation Size: No.10 / No.10A / No.11 / No.12 / No.14 / No.15 / No.15C / No.18 /No.20 / No.21 / No.22 / No.23 / No.24 / No.25 / No.36

OTHER MEDICAL CONSUMABLE



ALCOHOL PAD 70% Isopropyl alcohol



COTTON APPLICATOR 100% pure cotton





Changzhou Health Import and Export Company Ltd

Add: Sanhekou Street, Zhenglu Town, Tianning District, Changzhou, Jiangsu, China

Tel: +86-519-88670522 Fax: +86 519 88671958

E-Mail: info@tongdamedical.com

Website: www.healthtopmed.com









